The R&D Notebook

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1. Verity, I.; "Solvent Purity and the Individual Requirements of Different HPLC Detectors" LC-GC INT., Volume 8, Number 11, 652-655, 1995
2. Miyabe, K.; Takeuchi, S.; "Effect of Acetonitrile / Water Mobile Phase Composition on Adsorption Characteristics of Reversed Phase Liquid Chromatography" Analytical Chemistry, Volume 69, 2567-2574, 1997
3. Tanaka, M.; Yamazaki, H.; "Direct Determination of Pantoprazole Enantiomers in Human Serum by Reversed-Phase High-Performance Liquid Chromatography Using a Cellulose-Based Chiral Stationary Phase and Column Switching System as a Sample Cleanup Procedure" Analytical Chemistry, Volume 69, Number 9, 1513-1516, 1996
4. Kuangjing Shao, L.; Locke, D.C.; "Determination of Paclitaxel and Related Taxanes in Bulk Drug and Injectable Dosage Forms by Reversed Phase Liquid Chromatography" Analytical Chemistry, Volume 69, Number 11, 2008-2016, 1997
5. Corran, P.H.; "Reversed phase chromatography of proteins" in HPLC of Macromolecules. A practical approach. (ed. Oliver, R.W.A.) IRL Press, Oxford, 120-156, 1989
6. McCown, S.M.; Southern, D.; Morrison, B.E.; "Solvent properties and their effects on gradient elution highperformance liquid chromatography. III Experimental findings for water and acetonitrile" Journal of Chromatography, 352, 493-509, 1986
7. Darbouret, D.; Ishii, N.; Kanazawa, M.; Kano, I.; "Impact of feed water quality on the production of ultrapure water" The R&D Notebook, Millipore RD004 (2000)
8. Denoncourt, J.P.; Egozy, Y.; "Trace-level Analysis of High Purity Water Part III: High Performance Liquid Chromatography" Ultrapure Water January/February 1987
9. Anantharaman, V.; Parekh, B.; Hegde, R.; "Detection and Characterization of Organics in Semiconductor DI Water Processes" Ultrapure Water 11(3), 30-36, 1994
10. Darbouret, D.; Kano, I.; Youf, E.; Stewart, B.; "Optimizing Storage of purified water for Laboratory Applications" The R&D Notebook, Millipore RD001, 1998
11. Gabler, R.; Hegde, R.; Hughes, D.; "Degradation of High Purity Water on Storage" Journal of liquid Chromatography, 6 (13), 2565-2570
12. Achilli, G.; Cellerino, G.P.; Melzi d'Eril, G.; "Simultaneous determination of 27 phenols and herbicides in water by high-performance liquid chromatography with multielectrode electrochemical detection" Journal of Chromatography A, 697, 357-362, 1995
13. Coquart, V.; Hennion, M.C.; "Determination of Chlorotriazines in Aqueous Environmental Samples at the ng/l Level Using Preconcentration with a Cation Exchanger and On-Line High-Performance Liquid Chromatography" Journal of Chromatography, 585, 67-73, 1991
14. Somsen, G.W.; Hooyschuur, E.W.J.; Gooijer, C.; Brinkman, U.A.Th.; Velthorst, N.H.; Visser, T.; "Coupling of Reversed Phase Liquid Column Chromatography and Fourier Transform Infrared Spectrometry Using Postcolumn On-Line Extraction and Solvent Elimination" Analytical Chemistry, 68 (5), 746-752, 1996
15. Clark, K.; Retzik, M.; Darbouret, D. "Measuring TOC to maintain high-purity water" Ultrapure Water 14 (2), pp 21-24, (1997)

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